

## Defibtech Standalone Trainer

The Standalone AED Trainer can be used to instruct students in the various operating modes of Defibtech DDU-100 Series AEDs, and in the use of an AED as part of a standard rescue protocol. The DDU-100TR is similar to a functional DDU-100 series AED in form and operation, but all rubberized exterior surfaces are red, instead of black, to clearly distinguish the unit as a training device that cannot be used to defibrillate a real patient.

### Features

- Clearly identifiable as a training device
- Clear voice prompts identical to the real AED
- A remote control with 8 pre-programmed scenarios
- Pause button to allow for explanations
- Re-usable training pads

The Standalone AED Trainer is sold as part of the DCF-350T Standalone Training Package, which includes the following:

- DDU-100TR Trainer AED
- DBP-RC2 Training Battery (rechargeable)
- DTR-201 Training Battery Charger
- DDP-100TR Re-usable Adult Training Pads
- DTR-400 Remote Control
- Also available: DDP-200P-TR Re-usable Child/Infant Training Pads

### Training Remote Control

The remote control included can be used to control up to four standalone trainers or Lifeline AEDs with training battery packs installed. The remote control can be used to adjust speaker volume, for scenario selection, and for manual simulation of ventricular fibrillation, asystole, patient motion, normal sinus rhythm or incorrect pad placement.

